



## **XML ECG Decompressor** (available in English only)

Click on a blue button to select a file

---

This software is provided as-is, without any express or implied warranty. In no event shall Philips Medical Systems be held liable for any damages arising from the use of this software.

---

The XML ECG Decompressor application is used to convert any standard Philips ECG file in XML format (with compressed and base64 encoded ECG waveforms) into an ECG XML file with uncompressed, plain text, waveform sample values.

The XML ECG Decompressor application is compatible with any XML ECG file created on a PageWriter Trim or PageWriterTouch cardiograph, or with a Philips Intellivue Bedside Monitor and Information Center (with 12-lead capability), or with any XML file exported with the TraceMasterVue ECG Management System.

The XML ECG Decompressor application operates with Microsoft XML Core Services (MSXML 4.0 library) version 4.0 or higher. Click the button below for more information.

More  
Info

Click to read the program documentation as a Readme.txt file

More  
Info

Click to run the setup.exe installer

More  
Info

Click to read the program documentation for the Microsoft  
MSXML 4.0 Service Pack 2 program

---